


## PRODUCT SPECIFICATIONS

### Features:

- Auto-on and easy use
- Device pairing: pressing any of the devices, lights all three devices
- Disposable: discarded along with the infusion line
- Battery life: approximately 100 hours after activation

### Product Description

- Packaged in factory: three paired units per pack
- Pairing uses AES 128-bit encryption
- Operates at 915mHz on the ISM band for FCC Region 2 compliance
- Compatible with IV lines from 3.0mm to 4.2mm OD
- Model # IDA01



### General Data

- FDA Class I Exempt device
- Package dimension: 127mm x 230mm x 19mm
- Closed device dimensions: 61mm x 28mm x 17mm
- Operating temperature: -30 to 60 degrees C

### Materials Housing:

- Polypropylene Battery Chemistry: Lithium Manganese Dioxide (LiMnO<sub>2</sub>)