


PRODUCT SPECIFICATIONS

Features:

- Auto-on and easy use
- Device pairing: pressing any of the devices, lights all three devices
- Disposable: discarded along with the infusion line
- Battery life: at least 90 plus hours after activation

Product Description

- Packaged in factory: three paired units per pack
- Pairing uses AES 128-bit encryption
- Operates at 915mHz on the ISM band for FCC Region 2 compliance
- Compatible with IV lines from 3.0mm to 4.2mm OD
- Model # IDA01



General Data

- FDA Class I 510(k) Exempt device
- Package dimension: 127mm x 230mm x 19mm
- Closed device dimensions: 61mm x 28mm x 17mm
- Operating temperature: -30 to 60 degrees C

Material Housing and Battery:

- Housing Material: Polypropylene
- Battery Chemistry: Lithium Manganese Dioxide (LiMnO₂)
LiMnO₂ does not contain mercury, cadmium, or lead. Classified as non-hazardous waste according to Federal RCRA regulation (40 CFR 261) and are not regulated by DOT 49 CFR 173.185 (h) as hazardous materials. Product can be disposed by incineration as long as a controlled incinerator is utilized.

